

VISION-EASE LENS

WORLDWIDE

PRODUCT REFERENCE GUIDE

Polycarbonate

Advanta™ Spherical 65/75mm FSV 65mm/75mm Finished Single Vision

Equivalent Continua® Spherical SFSV Base Curve Selection

Sphere Power Up to -2.00 cyl	Base Curve	Curve Radius	Center Thickness	Edge Thickness
+4.00	6.33	83.728	4.9	1.2
+3.75	6.33	83.728	4.7	1.2
+3.50	6.33	83.728	4.5	1.2
+3.25	6.33	83.728	4.2	1.2
+3.00	6.33	83.728	4.0	1.2
+2.75	6.33	83.728	3.8	1.2
+2.50	6.33	83.728	3.5	1.2
+2.25	5.20	101.923	3.3	1.2
+2.00	5.20	101.923	3.0	1.2
+1.75	5.20	101.923	2.8	1.2
+1.50	5.20	101.923	2.6	1.2
+1.25	5.20	101.923	2.4	1.2
+1.00	5.20	101.923	2.1	1.2
+0.75	5.20	101.923	1.9	1.2
+0.50	5.20	101.923	1.9	1.4
+0.25	5.20	101.923	1.8	1.6
Plano	4.07	130.221	1.8	1.8
-0.25	4.07	130.221	1.6	1.9
-0.50	4.07	130.221	1.5	2.2
-0.75	4.07	130.221	1.5	2.5
-1.00	2.94	180.272	1.5	2.8
-1.25	2.94	180.272	1.5	3.1
-1.50	2.94	180.272	1.5	3.4
-1.75	2.94	180.272	1.4	3.7
-2.00	2.94	180.272	1.4	4.0
-2.25	2.94	180.272	1.4	4.3
-2.50	2.94	180.272	1.4	4.6
-2.75	2.94	180.272	1.4	5.0
-3.00	2.03	261.084	1.4	5.2
-3.25	2.03	261.084	1.3	5.5
-3.50	2.03	261.084	1.3	5.8
-3.75	2.03	261.084	1.3	6.1
-4.00	2.03	261.084	1.3	6.5
-4.25	2.03	261.084	1.3	6.8
-4.50	2.03	261.084	1.3	7.2
-4.75	2.03	261.084	1.3	7.5
-5.00	2.03	469.027	1.3	7.9
-5.25	1.13	469.027	1.3	8.0
-5.50	1.13	469.027	1.3	8.4
-5.75	1.13	469.027	1.3	8.8
-6.00	1.13	469.027	1.3	9.1
-6.25	1.13	469.027	1.3	9.5
-6.50	1.13	469.027	1.3	9.9
-6.75	1.13	469.027	1.3	10.2
-7.00	0.45	1177.778	1.3	10.4
-7.25	0.45	1177.778	1.3	10.8
-7.50	0.45	1177.778	1.3	11.1
-7.75	0.45	1177.778	1.3	11.5
-8.00	0.45	1177.778	1.3	11.9

Sphere	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00
65mm	4.00	6.25 6.33 – True Curve							
	3.75								
	3.50								
	3.25								
	3.00								
	2.75	5.25 5.20 – True Curve							
	2.50								
	2.25								
	2.00								
	1.75								
	1.50	4.00 4.07 – True Curve							
	1.25								
	1.00								
	0.75								
	0.50								
0.25	3.00 2.94 – True Curve								
0.00									
-0.25									
-0.50									
-0.75									
-1.00	2.00 2.03 – True Curve								
-1.25									
-1.50									
-1.75									
-2.00									
-2.25	1.00 1.13 – True Curve								
-2.50									
-2.75									
-3.00									
-3.25									
-3.50	0.50 0.45 – True Curve								
-3.75									
-4.00									
-4.25									
-4.50									
-4.75									
-5.00									
-5.25									
-5.50									
-5.75									
-6.00									
-6.25									
-6.50									
-6.75									
-7.00									
-7.25									
-7.50									
-7.75									
-8.00									
75mm	0.00	4.00 4.07 – True Curve							
	-0.25								
	-0.50								
	-0.75								
	-1.00								
	-1.25	3.00 2.94 – True Curve							
	-1.50								
	-1.75								
	-2.00								
	-2.25								
	-2.50	2.00 2.03 – True Curve							
	-2.75								
	-3.00								
	-3.25								
	-3.50								
-3.75	1.00 1.13 – True Curve								
-4.00									
-4.25									
-4.50									
-4.75									
-5.00	0.50 0.45 – True Curve								
-5.25									
-5.50									
-5.75									
-6.00									
-6.25									
-6.50									
-6.75									
-7.00									
-7.25									
-7.50									
-7.75									
-8.00									