


Product Reference Guide

Polycarbonate

VELA Poly FSV Super-Hydro AR
VELA Poly FSV AR
 65mm, 73mm Finished Single Vision

Sphere Power Up to -4.00 cyl	Base Curve	CT/ET
4.00	5.00	6.2/1.3
3.75	4.75	5.7/1.3
3.50	4.50	5.4/1.3
3.25	4.25	4.7/1.3
3.00	5.00	4.6/1.3
2.75	4.75	4.3/1.3
2.50	4.50	4.1/1.3
2.25	4.25	3.7/1.3
2.00	4.00	3.7/1.5
1.75	3.75	3.1/1.5
1.50	3.50	3.0/1.5
1.25	3.25	2.8/1.5
1.00	3.00	2.4/1.5
0.75	2.75	2.2/1.5
0.50	2.50	1.9/1.5
0.25	2.25	1.7/1.5
Plano	5.00	1.8
-0.25	5.00	1.7
-0.50	5.00	1.7
-0.75	5.00	1.7
-1.00	5.00	1.7
-1.25	5.00	1.7
-1.50	5.00	1.7
-1.75	5.00	1.7
-2.00	3.00	1.5
-2.25	3.00	1.5
-2.50	3.00	1.5
-2.75	3.00	1.5
-3.00	3.00	1.5
-3.25	3.00	1.5
-3.50	3.00	1.5
-3.75	3.00	1.5
-4.00	1.00	1.4
-4.25	1.00	1.4
-4.50	1.00	1.4
-4.75	1.00	1.4
-5.00	1.00	1.4
-5.25	1.00	1.4
-5.50	1.00	1.4
-5.75	1.00	1.4
-6.00	1.00	1.4

Base Curve Selection



Product Reference Guide

Polycarbonate

VELA Poly FSV

65mm, 73mm Finished Single Vision

Base Curve Selection

Sphere Power Up to -4.00 cyl	Base Curve	CT/ET
4.00	5.00	6.2/1.3
3.75	4.75	5.7/1.3
3.50	4.50	5.4/1.3
3.25	4.25	4.7/1.3
3.00	5.00	4.6/1.3
2.75	4.75	4.3/1.3
2.50	4.50	4.1/1.3
2.25	4.25	3.7/1.3
2.00	4.00	3.7/1.5
1.75	3.75	3.1/1.5
1.50	3.50	3.0/1.5
1.25	3.25	2.8/1.5
1.00	3.00	2.4/1.5
0.75	2.75	2.2/1.5
0.50	2.50	1.9/1.5
0.25	2.25	1.7/1.5
Plano	5.00	1.8
-0.25	5.00	1.7
-0.50	5.00	1.7
-0.75	5.00	1.7
-1.00	5.00	1.7
-1.25	5.00	1.7
-1.50	5.00	1.7
-1.75	5.00	1.7
-2.00	3.00	1.5
-2.25	3.00	1.5
-2.50	3.00	1.5
-2.75	3.00	1.5
-3.00	3.00	1.5
-3.25	3.00	1.5
-3.50	3.00	1.5
-3.75	3.00	1.5
-4.00	1.00	1.4
-4.25	1.00	1.4
-4.50	1.00	1.4
-4.75	1.00	1.4
-5.00	1.00	1.4
-5.25	1.00	1.4
-5.50	1.00	1.4
-5.75	1.00	1.4
-6.00	1.00	1.4

