


Product Reference Guide



Polycarbonate

VISION EASE Polarized Tints
Semi-Finished Single Vision
 75mm Polarized Single Vision (Gray 2)

Base Curve	True Curve (1.530)	Front Radius	50mm SAG*	Back Curve	Back Radius	CT (mm) Nominal	CT (mm) Minimum	ET (mm) Nominal
2.00	2.03	261.08	1.20	-6.25	-84.80	8.0	7.5	14.0
4.00	4.07	130.22	2.42	-6.25	-84.80	8.8	8.3	12.0
6.25	6.33	83.73	3.82	-6.25	-84.80	8.6	8.1	8.5
8.50	8.60	61.63	5.30	-7.50	-70.67	8.9	8.4	7.0

*Utilized for power calculation only, not for thickness calculation. May involve limitations due to machine configuration.

Base Curve Selection

